



QUARTET
TECHNOLOGY, INC.

Radio Remote(s)
(Installation Manual)

87 Progress Avenue, Tyngsboro, MA 01879, USA
Phone (978) 649-4ECU Fax (978) 649-8363

<http://www.qtiusa.com>

Trademarks, Version, Printing, and Copyright

Trademarks

Simplicity™ All-In-One™, Simplicity™ Voice, Simplicity™ Switch, Smart Switch™, VMK™, and QSP™ are registered trademarks of Quartet Technology, Inc.

Version

August, 2001
Version 1.0
North American Edition

Printing

Printed in the USA.

Quartet Technology, Inc. (QTI) is concerned about the environment. To reduce waste and complete the recycling circle, we printed this manual and cover on stock that is recyclable. QTI has made every effort to look at environmental implications when deciding on packaging.

Copyright

Copyright© 2001 Quartet Technology, Inc.

Radio Remote(s) Installation Manual P/N: 5093

This manual is copyrighted and all rights are reserved. This manual may not, in whole or in part, be copied, photocopied, translated, or reduced to any electronic medium or machine-readable form without prior consent, in writing, from Quartet Technology, Inc.

The information in this manual is subject to change without notice and Quartet Technology, Inc. assumes no responsibility for any errors that may appear in this document.

Warranty Information

Warranty

In connection with Products sold and Software licensed by Quartet to the Buyer hereunder, Quartet makes the following warranties, and no others:

(1) Quartet warrants that the Products (excluding, where applicable, any batteries, microphones, therein contained, or any portion thereof comprised of Software, which is specifically warranted below) will be free of defects in materials and workmanship for a period of one year from the date of delivery to the Buyer. Quartet's sole responsibility under the aforesaid warranty shall be, at its election, the repair or replacement of defective materials and/or workmanship during the aforesaid period, provided that the Buyer has promptly reported such defect to Quartet and Quartet has found, upon inspection, that such Products are defective. The Buyer must obtain shipping instructions from Quartet prior to returning any Products under warranty. It is the Buyer's responsibility to return the defective Products to Quartet at the Buyer's expense. Ordinary transportation charges back to the Buyer shall be the responsibility of the Buyer, exclusive of duties, tariffs and transportation insurance. All replaced hardware or parts shall become Quartet's property.

(2) Quartet warrants that all unmodified Software will substantially conform, for a period of one year from the date of delivery to the Buyer, to Quartet's published specifications prevailing at the time of shipment. Receipt by Quartet of a written claim under this warranty and the return of the Software must occur within said period. Quartet does not warrant that the Software will be error free or that all errors will be remedied. Quartet's entire liability under this warranty shall be for Quartet to make reasonable efforts to remedy, in a manner deemed suitable to Quartet, any nonconformance reported in writing during the warranty period.

(3) EXCEPT FOR THE WARRANTIES SET FORTH IN THIS SUBSECTION, QUARTET HEREBY DISCLAIMS, TO THE MAXIMUM EXTENT ALLOWED BY LAW, ALL OTHER EXPRESS OR IMPLIED WARRANTIES IN RESPECT OF THE PRODUCTS OR THE SOFTWARE, INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTIES OF MERCHANTABILITY AND/OR FITNESS FOR A PARTICULAR PURPOSE. THE EXPRESS WARRANTIES STATED ABOVE ARE IN LIEU OF ALL OBLIGATIONS AND LIABILITIES OF QUARTET FOR DAMAGES, INCLUDING, WITHOUT LIMITATION, INCIDENTAL CONSEQUENTIAL (INCLUDING, WITHOUT LIMITATION, LOSS OF PROFITS OR LOSS OR DESTRUCTION OF DATA), INDIRECT, SPECIAL OR EXEMPLARY DAMAGES THAT ARISE IN ANY WAY OUT OF ANY DEFECTS IN THE PRODUCTS OR SOFTWARE OR IN THEIR MATERIALS OR WORKMANSHIP

Manual Description

Audience

This manual is for users of any of the following Quartet Radio Remote packages:

QTI P/N:	Description
9219	Remote Switch Package
9221	Remote Voice Package
9231	Remote AIO Package
9250	Voice Back Radio Package
9261	Remote Switch Package w/Battery Adapter
9262	Remote Voice Package w/Battery Adapter
9263	Remote AIO Package w/Battery Adapter

Purpose

The purpose of this manual is to provide:

- Identification of components
- Proper connections
- Suggested mounting

Intended Use

The intended function of the Quartet radio remote packages is to provide one way and/or two way communication with Quartet's Simplicity line of ECU's using voice, switch, or voice and switch activation.

(This page left blank intentionally)

Table of Contents

TRADEMARKS, VERSION, PRINTING, AND COPYRIGHT	2
WARRANTY INFORMATION	3
MANUAL DESCRIPTION	4
OVERVIEW	8
<i>Introduction</i>	8
<i>Features</i>	8
<i>Objectives</i>	8
<i>Caution</i>	8
REMOTE SWITCH PACKAGE (P/N: 9219)	9
<i>Transmitter View</i>	9
<i>Transmitter Components</i>	9
<i>Receiver View</i>	10
<i>Receiver Components</i>	10
REMOTE VOICE PACKAGE (P/N: 9221)	11
<i>Transmitter View</i>	11
<i>Transmitter Components</i>	11
<i>Receiver View</i>	12
<i>Receiver Components</i>	12
REMOTE AIO PACKAGE (P/N: 9231)	13
<i>Transmitter View</i>	13
<i>Transmitter Components</i>	13
<i>Receiver View</i>	14
<i>Receiver Components</i>	14
VOICE BACK RADIO PACKAGE (P/N: 9250)	15
<i>Transmitter View</i>	15
<i>Transmitter Components</i>	15
<i>Receiver View</i>	16
<i>Receiver Components</i>	16
<i>Receiver Power Connection</i>	16
<i>Receiver Controls</i>	17
REMOTE SWITCH PACKAGE W/ BATTERY ADAPTER (P/N: 9261)	18
<i>Transmitter View</i>	18
<i>Transmitter Components</i>	19
<i>Receiver View</i>	20
<i>Receiver Components</i>	20
REMOTE VOICE PACKAGE W/ BATTERY ADAPTER (P/N: 9262)	21
<i>Transmitter View</i>	21
<i>Transmitter Components</i>	21
<i>Receiver View</i>	22
<i>Receiver Components</i>	22
REMOTE AIO PACKAGE W/ BATTERY ADAPTER (P/N: 9263)	23
<i>Transmitter View</i>	23
<i>Transmitter Components</i>	24
<i>Receiver View</i>	25
<i>Receiver Components</i>	25

Table of Contents, *continued*

APPENDIX A	26
<i>Connecting Switches to the Remote Encoder</i>	27
APPENDIX B	29
<i>Wheelchair Pouch Component Placement</i>	30
APPENDIX C	32
<i>Wheelchair Mounting</i>	33
APPENDIX D	34
<i>Battery Charger</i>	35
APPENDIX E	36
<i>Transmitter “On/Off” Switch</i>	37
<i>Transmitter View</i>	37
<i>Transmitter Components</i>	37
APPENDIX F	38
<i>Cleaning and Maintenance</i>	39

Overview

Introduction

Congratulations on your purchase of a Quartet radio remote package. By adding a radio remote package to your Simplicity ECU, you have enhanced its capabilities and gained increased independence by being free to roam within your environment.

Features

Features of the radio remote package are:

- 150ft operating range
 - UHF transmission band
 - Crystal controlled frequency (no adjustments)
 - True diversity reception
-

Objectives

In this manual you will learn about:

- Identifying components for each package
 - Connecting the components
 - Placement of components in the wheelchair pouch
 - Mounting of components on the wheelchair
 - Receiver placement tips
 - Battery Charger
 - Accessories
 - Cleaning and maintenance
-

Caution

Electrical shock can result from the removal of the receiver cover. Refer servicing to qualified service personnel. There are no user-serviceable parts inside. Do not expose to rain or moisture.


The circuits inside the transmitter and receiver have been precisely adjusted for optimum performance and compliance with federal regulations. Do not attempt to open the receiver or transmitter. To do so will void the warranty, and may cause improper operation.

Remote Switch Package (P/N: 9219)

Introduction

The **Remote Switch Package** (P/N: 9219) is used with Simplicity™ Switch ECU's. It provides remote switch capability (single and/or dual operation) to the ECU as well as remote microphone capability for use with the telephone or computer.

Transmitter View



Transmitter Components

Letter	Component	QTI P/N
A	Microphone	9245
B	Remote Switch Encoder	9240
C	Transmitter Cable	9117
D*	Switch Cable	9256
E*	QSP Switch	9260
F	Battery Pack	9080
G**	UHF Transmitter (single frequency)	9217

* Optional components

** The power switch **must** be left in the “on” position (*See Appendix E*)


Note

The transmitter includes a permanently attached flexible antenna. For best results, allow the antenna to hang freely and full length from the transmitter.

Do not attempt to remove, replace, or change the length of the transmitting antenna.

Remote Switch Package (P/N: 9219) (cont.)

Receiver View



Receiver Components

Letter	Component	QTI P/N
A	UHF Receiver AC Power Adapter	4312
B	Receiver Audio Cable	4049
C	UHF Receiver (single frequency)	9235

Note

The antennas on the receiver normally should be positioned in the shape of a “V” (45° from vertical) for best reception.

Note

The receiver has been carefully calibrated at the factory for optimum performance. The knob on the front of the receiver has been glued in place to prevent accidental movement.


Do not attempt to remove the adhesive or change the setting of this control.

Remote Voice Package (P/N: 9221)

Introduction

The **Remote Voice Package** (P/N: 9221) is used with Simplicity™ Voice ECU's. It provides remote voice capability to the ECU including audio for use in telephony and computers.

Transmitter View



Transmitter Components

Letter	Component	QTI P/N
A	Microphone	9184
B	Female to female adapter	4047
C	Transmitter Cable	9117
D**	UHF Transmitter (single frequency)	9217
E	Battery Pack	9080

** The power switch **must** be left in the “on” position (*See Appendix E*)


Note

The transmitter includes a permanently attached flexible antenna. For best results, allow the antenna to hang freely and full length from the transmitter.

Do not attempt to remove, replace, or change the length of the transmitting antenna.

Remote Voice Package (P/N: 9221) (cont.)

Receiver View



Receiver Components

Letter	Component	QTI P/N
A	UHF Receiver AC Power Adapter	4312
B	Receiver Audio Cable	4049
C	UHF Receiver (single frequency)	9235

Note

The antennas on the receiver normally should be positioned in the shape of a “V” (45° from vertical) for best reception.

Note

The receiver has been carefully calibrated at the factory for optimum performance. The knob on the front of the receiver has been glued in place to prevent accidental movement.


Do not attempt to remove the adhesive or change the setting of this control.

Remote AIO Package (P/N: 9231)

Introduction

The **Remote AIO Package** (P/N: 9231) is used with Simplicity™ AIO ECU's. It provides remote switch capability (single and/or dual operation) to the ECU as well as remote microphone capability for use with the telephone or computer.

Transmitter View



Transmitter Components

Letter	Component	QTI P/N
A	Microphone	9245
B	Remote Switch Encoder	9240
C	Transmitter Cable	9117
D*	QSP Switch Cable (1')	9256
E*	QSP Switch	9260
F	Battery Pack	9080
G**	UHF Transmitter (single frequency)	9217

* Optional components

** The power switch **must** be left in the “on” position (*See Appendix E*)


Note

The transmitter includes a permanently attached flexible antenna. For best results, allow the antenna to hang freely and full length from the transmitter.

Do not attempt to remove, replace, or change the length of the transmitting antenna.

Remote AIO Package (P/N: 9231) (cont.)

Receiver View



Receiver Components

Letter	Component	QTI P/N
A	UHF Receiver AC Power Adapter	4312
B	Receiver Audio Cable	4049
C	UHF Receiver (single frequency)	9235

Note

The antennas on the receiver normally should be positioned in the shape of a “V” (45° from vertical) for best reception.

Note

The receiver has been carefully calibrated at the factory for optimum performance. The knob on the front of the receiver has been glued in place to prevent accidental movement.


Do not attempt to remove the adhesive or change the setting of this control.

Voice Back Radio Package (P/N: 9250)

Introduction

The **Voice Back Radio Package** (P/N: 9250) can be used with all Simplicity™ ECU's. It provides remote listening capability of the ECU.

Transmitter View



Transmitter Components

Letter	Component	QTI P/N
A**	UHF Transmitter (100-Channel) w/Power Cord	9222
B	Transmitter Cable	9124
C	UHF Transmitter AC Power Adapter	8087

** The power switch **must** be left in the “on” position (*See Appendix E*)


Note

The transmitter includes a permanently attached flexible antenna. For best results, allow the antenna to extend freely and full length from the transmitter.

Do not attempt to remove, replace, or change the length of the transmitting antenna.

Voice Back Radio Package (P/N: 9250) (cont.)

Receiver View




Receiver Components

Letter	Component	QTI P/N
A	Mini-Cube Speaker	8086
B	UHF Receiver (100-Channel)	4295

Note

Use the “Monitor Level” knob to adjust the listening volume.

Receiver Power Connection



Note


Depending on the particular radio remote package used, the connector labeled “A” in the above illustration will connect either to the transmitter or remote encoder.

Remote Switch Package w/Battery Adapter (P/N:9261)

Introduction

The *Remote Switch Package w/Battery Adapter* (P/N: 9261) is used with Simplicity™ Switch ECU's. It provides remote switch capability (single and/or dual operation) to the ECU as well as remote microphone capability for use with the telephone or computer.

Transmitter View



Remote Switch Package w/Battery Adapter (P/N:9261) (cont.)

Transmitter Components

Letter	Component	QTI P/N
A	Battery Adapter	9248
B	Battery Adapter +9VDC Output Power Cable	9254
C	Battery Adapter Wheelchair Power Cable	9252
D**	UHF Transmitter (single frequency)	9217
E	Transmitter Cable	9117
F	Microphone	9245
G*	QSP Switch	9260
H*	QSP Switch Cable (1')	9256
I	Remote Switch Encoder	9240

* Optional components

** The power switch **must** be left in the “on” position (*See Appendix E*)


Note

The transmitter includes a permanently attached flexible antenna. For best results, allow the antenna to hang freely and full length from the bottom of the transmitter.

Do not attempt to remove, replace, or change the length of the transmitting antenna.

Remote Switch Package w/Battery Adapter (P/N:9261) (cont.)

Receiver View



Receiver Components

Letter	Component	QTI P/N
A	UHF Receiver AC Power Adapter	4312
B	Receiver Audio Cable	4049
C	UHF Receiver (single frequency)	9235

Note

The antennas on the receiver normally should be positioned in the shape of a “V” (45° from vertical) for best reception.

Note

The receiver has been carefully calibrated at the factory for optimum performance. The knob on the front of the receiver has been glued in place to prevent accidental movement.


Do not attempt to remove the adhesive or change the setting of this control.

Remote Voice Package w/Battery Adapter (P/N:9262)

Introduction

The **Remote Voice Package w/Battery Adapter** (P/N: 9262) is used with Simplicity™ Switch ECU's. It provides remote voice capability to the ECU including audio for use with the telephone or computer.

Transmitter View



Transmitter Components

Letter	Component	QTI P/N
A	Battery Adapter	9248
B	Battery Adapter +9VDC Output Power Cable	9254
C	Battery Adapter Wheelchair Power Cable	9252
D**	UHF Transmitter (single frequency)	9217
E	Transmitter Cable	9117
F	Female-to-Female Coupler	4047
G	Microphone	9184

** The power switch **must** be left in the “on” position (*See Appendix E*)


Note

The transmitter includes a permanently attached flexible antenna. For best results, allow the antenna to hang freely and full length from the bottom of the transmitter.

Do not attempt to remove, replace, or change the length of the transmitting antenna.

Remote Voice Package w/Battery Adapter (P/N:9262) (cont.)

Receiver View



Receiver Components

Letter	Component	QTI P/N
A	UHF Receiver AC Power Adapter	4312
B	Receiver Audio Cable	4049
C	UHF Receiver (single frequency)	9235

Note

The antennas on the receiver normally should be positioned in the shape of a “V” (45° from vertical) for best reception.

Note

The receiver has been carefully calibrated at the factory for optimum performance. The knob on the front of the receiver has been glued in place to prevent accidental movement.


Do not attempt to remove the adhesive or change the setting of this control.

Remote AIO Package w/Battery Adapter (P/N:9263)

Introduction

The *Remote AIO Package w/Battery Adapter* (P/N: 9263) is used with Simplicity™ Switch ECU's. It provides remote switch capability (single and/or dual operation) to the ECU as well as remote microphone capability for use in telephony and computers.

Transmitter View



Remote AIO Package w/Battery Adapter (P/N:9263) (cont.)

Transmitter Components

Letter	Component	QTI P/N
A	Battery Adapter	9248
B	Battery Adapter +9VDC Output Power Cable	9254
C	Battery Adapter Wheelchair Power Cable	9252
D**	UHF Transmitter (single frequency)	9217
E	Transmitter Cable	9240
F	Microphone	9184
G*	QSP Switch	9260
H*	QSP Switch Cable (1')	9256
I	Remote AIO Encoder	9239

* Optional components

** The power switch **must** be left in the “on” position (*See Appendix E*)


Note

The transmitter includes a permanently attached flexible antenna. For best results, allow the antenna to hang freely and full length from the bottom of the transmitter.

Do not attempt to remove, replace, or change the length of the transmitting antenna.

Remote AIO Package w/Battery Adapter (P/N:9263) (cont.)

Receiver View



Receiver Components

Letter	Component	QTI P/N
A	UHF Receiver AC Power Adapter	4312
B	Receiver Audio Cable	4049
C	UHF Receiver (single frequency)	9235

Note

The antennas on the receiver normally should be positioned in the shape of a “V” (45° from vertical) for best reception.

Note

The receiver has been carefully calibrated at the factory for optimum performance. The knob on the front of the receiver has been glued in place to prevent accidental movement.

Do not attempt to remove the adhesive or change the setting of this control.

Appendix A

Connecting Switches to the Remote Encoder

Introduction


Both the *Remote Switch Package* and the *Remote AIO Package* contain a device called the remote encoder. As its name implies, this device encodes both the voice and switch signals into one audio channel that is transmitted to the ECU.

When connecting any switch (single or dual) to the encoder, you **must** use a stereo plug.

Stereo Plugs

Before connecting a switch (single or dual) to the remote encoder, it is important to understand what a stereo plug is and the terminology associated with it.

Shown below is a cross-section of a typical stereo plug.



Note

All references in this manual that refer to tip, ring, and sleeve assume the connections as shown above.

Note

Almost all adaptive devices require a connection between tip and sleeve, ring and sleeve, or both.

Connecting Switches to the Remote Encoder (cont.)

Switch Combinations

The following table illustrates the various single and dual switch combinations between the ECU and the remote radio packages.

<i>ECU</i>	<i>Remote Encoder</i>
Dual Switch (Stereo Plug)	Dual Switch (Stereo Plug)
Dual Switch (Stereo Plug)	Single Switch (Stereo Plug)
Single Switch (Stereo Plug)	Dual Switch (Stereo Plug)
Single Switch (Mono Plug)	Single Switch (Stereo Plug)

Note

In all cases, when attaching any type of switch to the remote encoder, you ***must*** use a stereo plug.

Note

The Simplicity™ line of ECU's will automatically detect the type of plug connected to the ECU.

Note

Never use adapters or extension cables when connecting switches or microphones to the remote encoder.


Appendix B

Wheelchair Pouch Component Placement


Introduction

All radio remote packages come with a wheelchair pouch (P/N: 8085). The following illustrations show how the various components from the different radio remote packages are organized in the pouch.

Remote Switch Package w/ Battery Adapter (P/N: 9261)




Remote Voice Package w/ Battery Adapter (P/N: 9262)




Wheelchair Pouch Component Placement (cont.)

Remote AIO Package w/ Battery Adapter (P/N: 9263)



Remote Packages w/ Battery Adapter and Voice Back Radio




Appendix C

Wheelchair Mounting

Introduction


The following illustrations show how the various components from the different radio remote packages are typically mounted to a wheelchair.

Wheelchair Pouch Mounting



LA8183

Microphone, QSP Straw, and/or Mini-Cube Speaker Mounting



LA8182

Appendix D


Battery Charger

Introduction

The **Remote Switch Package** (P/N: 9219), **Remote Voice Package** (P/N: 9221), and the **Remote AIO Package** (P/N: 9231) all come with a battery charger and two rechargeable Ni-Cad batteries.

The batteries provide approximately 12 hours of use after a full charge. While one battery is in use, the other battery should be charging. The charger takes approximately 14 hours to charge a fully discharged battery.

Charger View



Charger Components

Letter	Component	QTI P/N
A	Ni-Cad Charger	9178
B	Ni-Cad Battery Assembly	9080

Note

The charger is automatic. The battery can be left connected to the charger indefinitely.

Note

It is best to alternate the batteries everyday by using one battery while the other is charging.


Appendix E

Transmitter On/Off Switch

Introduction

The transmitter power switch has three positions: *Off*, *Standby*, and *On*. The transmitter **must** be left in the “On” position at all times in order to function properly.

Transmitter View



Transmitter Components

Letter	Description
A	Power Switch
B	Power Indicator

Note

The power switch must be left in the “ON” position.

Note

For single channel transmitters (P/N: 9217), the power indicator LED will remain lit when the power switch is in the “**ON**” position.

Note

For multi-channel transmitters (P/N: 9222), the power indicator LED will flash on briefly when the power is switched to the “**ON**” position.

Appendix F

Cleaning and Maintenance

Cleaning

You can clean any of the radio components using a slightly damp cloth, when necessary. Use water only.

Caution

Do not get any liquid or foreign material in the components.

Maintenance

The radio packages require no special cleaning or daily maintenance.

Warning

You should check the radio package connections monthly for safety and integrity. Contact your authorized Quartet distributor for necessary repairs or replacements.

(This page left blank intentionally)